



Handheld Pulse Oximeter



Model: A360

Standard Accessory



Adult finger
SpO₂ sensor



USB line

Optional Accessory



Neonatal finger
SpO₂ sensor



Child finger
SpO₂ sensor



Temperature
sensor



Feature

- Silicon rubber shell protection and stable bracket for table usage.
- Using DB9 type connector compatible with Nellcor Spo2 Sensor.
- Support the oximeter probe for adult, neonate and infant.
- Big size 2.8 inch color TFT LCD display for SPO₂/PR/Pulse bargraph/ plethysmogram.
- Inner Flash memory can store testing result up to 300 hours.
- USB interface support upload the data to computer and review the history data with software in PC.
- Visual and sound alarm function.
- Adjust the parameters in friendly menu.
- Low Battery voltage indicator.
- Automatically switch off within 3 minutes when no signal.
- 15-hours continuous operation with a new, fully charged battery.

Specifications

Power

4 × AAA 1.5V Alkaline Battery or
rechargeable 1000mAh Li battery

Hemoglobin saturation display: 35%-100%

Pulse rate Display: 30-250BPM

Power Consumption:

smaller than 80mA (Normal)

Warranty: main unit (1 Year)

SpO₂ sensor (6 months)

Dimension (L×W×H): 145mm×83mm×33mm

Color box Dimension (L×W×H): 185mm×126mm×89mm

Resolution:

Hemoglobin Saturation(SpO₂): 1%
Pulse rate : 1BPM

Measurement Accuracy:

Hemoglobin Saturation(SpO₂):
2%(70%-100 %);unspecified(< 70%)

Pulse rate: 1BPM

Low perfusion Condition: 0.3%

Net Weight: 240 g

Gross Weight: 510 g