









# 8F Oxygen Concentrator

## 8F-5AW

-  95.5~87% Oxygen Concentration
-  U.S. Imported Sieve
-  Ultra Silent
-  High Reverse Display LCD Screen
-  Intelligent Surveillance
-  Infrared Remote Control System



Flow Range: 0 ~ 5L/min

Oxygen Concentration: 95.5 ~ 87%

Maximum Outlet Pressure: 40~70 kPa

Pressure Relief Mechanism Operational at: 250kPa ± 25kPa (36.25psi ± 3.63psi)

Altitude: Up to 1828 meters (6000 ft) above sea level without degradation of concentration levels.

From 1828 meters (6000 ft) to 4000meters (13129 ft) below 90% efficiency.

Sound Level: Front side 48dB(A)(average), Whole unit 52dB(A)

Maximum nebulizing rate (apply in the oxygen concentrators having the function of nebulizing): ≥0.15 mL/min

Power Supply:

AC220V ± 10%

50Hz ± 2%

Input Power: 400VA

Weight: 15.5kg (34.1 lbs)

Dimension: W39 cm x D 24.5cm x H50 cm

Minimum Operating Time: 30minutes

Oxygen outlet adapter is made of chromium plated copper, not become rusty

### JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Sales Center Address: Huanyuan East Road No.1, Xuzhuang Software Park, Nanjing, Jiangsu Province, P.R. China

Tel: 86-25-8713 6530 Fax: 86-25-8713 6530

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn

www.yuwell.com



/globalyuwell