




# Desktop Fetal Doppler



Model: A100S



Curve mode



Replay mode



## Feature

- OLED display.
- 6 Modes.
- 36 Hours FHR data storage and reload ability.
- High Sensitive 2.5 MHz Doppler probe, accurate FHR detection.
- RTC Timer inside can display date and time.
- Built-in rechargeable battery for 12 hours continuous work.
- Medical-use switch power.
- Support multinational language.

## Specifications

### FHR

Ultrasound Frequency: 2.5MHz  
Intensity: 2.5mW/cm or less  
FHR Range: 30-240 bpm  
FHR Accuracy:  $\pm 1$  bpm

### Function

Memory and Replay function  
Alarm on/off function  
Volume +/- function  
Menu Set function

### Display

OLED displays FHR number and curve,  
working mode, menu setting, date & time.  
Battery capacity display  
Power off/on display  
Alarm on/off display  
Probe off display

### Power Supply

Input: AC 100 - 240V (50/60 Hz)

### Standard Accessory

|                                     |      |
|-------------------------------------|------|
| 2.5MHz Doppler Probe                | 1 pc |
| Rechargeable Li-ion Battery 2200mAh | 1 pc |
| Ultrasound Gel                      | 1 pc |
| Operation Manual                    | 1 pc |
| Power Cord                          | 1pc  |
| GEL/Probe Bracket                   | 1pc  |

### Optional Accessory

|                      |     |
|----------------------|-----|
| 2.0MHz Doppler Probe | 1pc |
| Ultrasound Gel       | 1pc |

### Dimension:

280mm(L) x 148mm(W) x 200mm(H)

### Net Weight:

1.5 kg

### Outer box Dimension:

330mm (L) x 250mm (W) x 260mm (H)

### Gross weight:

3.2kg

### Warranty:

2 years (probe excluded)