




SKB2A12

SPINE BOARD




SKB2A12

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 1840×450×55 mm
- Load Bearing 159kg



Technical configuration

- Belt 3pcs

Material

This product is blow molded by high density plastic polyethylene, strong and durable, not easy to aging.




Function

Compatible with X-ray, MRI, CT transparent, convenient for examination of patient on board, reduce the possible hurt to unload the patient from the board.



Handle

The handle whole is evenly opened around, which can be used for many people to lift, carry and lift at the same time.



Float

Hard structure, easy to continue CPR and cardiac compression rescue in the process of transport, large displacement, the overall volume of 0.04m³, can float an adult in normal temperature water, greatly reducing the difficulty of water rescue personnel, convenient for faster rescue of the injured.



Bandages

This product is equipped with three bandages as standard. The optional kit includes neck bracket, head immobilizer and spine fix plate.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue